



**Pumani**  
bubble CPAP



MAX PRESSURE  
34.5 kPa

DO NOT CLOSE



**Pumani**  
bubble CPAP

Product Catalogue © 2018





# Table of Contents

## Pumani bubbleCPAP..... Page 3-5

What is included in the package?

What is included in the “Starter Kit”?

## Pumani Accessory Items..... Page 6-11

What are the various patient accessories?

How to order patient accessories and spare items?

## Product Recommendations..... Page 12-14

Nasal prong and hat size guide

Recommended accessory order quantities

# Pumani bubble CPAP



# What is included in the package?

The Pumani bubble CPAP package has:

- One Pumani bubble CPAP device
- Pumani 'Starter kit' (see next page)  
(Item ID 3-597)



# What is included in the 'Starter Kit'?



## Patient Tubing

Quantity: 3



## Bottle

Quantity: 1



## Fish Elbow

Quantity: 6



## Hat Clips

Quantity: 6



## Hats (Small, Medium, Large)

Quantity: 3

Small - 1      Medium - 1      Large - 1



## Variety Nasal Prong Pack

Quantity: 6 prongs, 1 each

Size 0 - 1      Size 3 - 1  
Size 1 - 1      Size 4 - 1  
Size 2 - 1      Size 5 - 1



## Tube Holder

Quantity: 3



## Power Cord, UK Style

Quantity: 1



## User Manual

Quantity: 1

## To Re-Order:

Re-order Pumani bubble CPAP device using Item ID, 3-597. Contact us:

[sales@pumani.com](mailto:sales@pumani.com)

[www.pumani.com](http://www.pumani.com)

+1 (628) 236 7335

### Please Note:

'Starter Kit' cannot be re-ordered. Please order "Patient Accessories" or "Spare Items" from page 6-10.

# Pumani Accessory Items



# Pumani Patient Accessories

There are 6 unique patient accessories as shown below. Each of these items can be ordered individually.



## To Order Patient Accessories:

To order these patient accessories,

Find out the "Item ID" for the particular patient accessory and contact us:

✉ [sales@pumani.com](mailto:sales@pumani.com)

🌐 [www.pumani.com](http://www.pumani.com)

☎ +1 (628) 236 7335

# Patient Accessories - Nasal Prong

Nasal Prongs can be ordered in single size packs or variety size packs.

## To Order Patient Accessories:

To order these patient accessories,

Find out the "Item ID" for the particular patient accessory and contact us:

✉ [sales@pumani.com](mailto:sales@pumani.com)

🌐 [www.pumani.com](http://www.pumani.com)

☎ +1 (628) 236 7335

Refer to page 13 for Nasal Prong Size Guide



000 Size 000




Nasal Prong  
(Size 000)

Pack of 6

Item ID: 3-548

0 Size 0



Nasal Prong  
(Size 0)

Pack of 6

Item ID: 3-520

1 Size 1



Nasal Prong  
(Size 1)

Pack of 6

Item ID: 3-521

2 Size 2



Nasal Prong  
(Size 2)

Pack of 6

Item ID: 3-522

3 Size 3



Nasal Prong  
(Size 3)

Pack of 6

Item ID: 3-523

4 Size 4



Nasal Prong  
(Size 4)

Pack of 6

Item ID: 3-524

5 Size 5



Nasal Prong  
(Size 5)

Pack of 6

Item ID: 3-526

0-5 Size 0-5  
(variety)



Nasal Prong  
Variety (Size 0-5)

Pack of 6 (1 each)

Item ID: 3-527



# Patient Accessories - Hat

## To Order Patient Accessories:

To order these patient accessories,

Find out the "Item ID" for the particular patient accessory and contact us:

✉ [sales@pumani.com](mailto:sales@pumani.com)

🌐 [www.pumani.com](http://www.pumani.com)

☎ +1 (628) 236 7335

Refer Page 13 for Hat Size Guide



Pumani Hat, Size - Small

Pack of 10

Item ID: 3-533



Pumani Hat, Size - Medium

Pack of 10

Item ID: 3-534



Pumani Hat, Size - Large

Pack of 10

Item ID: 3-535

# Patient Accessories Continued.,



## Hat Clip Pack

Pack of 10

Item ID: 3-536



## Tube Holder Pack

Pack of 10

Item ID: 3-552



## Fish Elbow Pack

Pack of 10

Item ID: 3-530



## Patient Tubing

1 set

Item ID: 3-551

## To Order Patient Accessories:

To order these patient accessories,

Find out the "Item ID" for the particular patient accessory and contact us:

✉ [sales@pumani.com](mailto:sales@pumani.com)

🌐 [www.pumani.com](http://www.pumani.com)

☎ +1 (628) 236 7335

# Pumani Spare Items



## Bottle

1 piece

Item ID: 3-549



## Filter

1 piece

Item ID: 3-541



## Flow Meter- Oxygen

1 piece

Item ID: 3-543



## Flow Meter-Total

1 piece

Item ID: 3-544

## To Order Patient Accessories:

To order these patient accessories,

Find out the "Item ID" for the particular patient accessory and contact us:

✉ [sales@pumani.com](mailto:sales@pumani.com)

🌐 [www.pumani.com](http://www.pumani.com)

☎ +1 (628) 236 7335

# Product Recommendations



# Nasal Prong Size and Hat Size Recommendations



### Pumani Nasal Prong Size Selection Chart

Patient Weight Range	Nasal Prong Size
Less than 1,000 grams	000 or 0
1,000 grams to 1,250 grams	1
1,250 grams to 2,000 grams	2
2,000 grams to 3,000 grams	3
3,000 grams to 4,000 grams	4
Over 4,000 grams	5

### Pumani Hat Size Selection Chart

Patient Weight Range	Hat Size
Less than 1,500 grams	Small
1,500 grams to 3,000 grams	Medium
Over 3,000 grams	Large

# Suggested Accessory Order Quantities

Accessory Items	Item #	10 Patients	25 Patients	Requested Quantity
Nasal Prong Pack, Size 000 (Pack of 6)*	3-548	1	1	
Nasal Prong Pack, Size 0 (Pack of 6)*	3-520	1	1	
Nasal Prong Pack, Size 1 (Pack of 6)*	3-521	1	1	
Nasal Prong Pack, Size 2 (Pack of 6)*	3-522	1	1	
Nasal Prong Pack, Size 3 (Pack of 6)*	3-523	1	1	
Nasal Prong Pack, Size 4 (Pack of 6)*	3-524		1	
Nasal Prong Pack, Size 5 (Pack of 6)*	3-526			
Nasal Prong Pack, Variety (Pack of 6)*	3-527	1	1	
Small Hat (Pack of 10)**	3-533	1	2	
Medium Hat (Pack of 10)**	3-534	1	2	
Large Hat (Pack of 10)**	3-535	1	1	
Hat Clip (Pack of 10)	3-536	1	1	
Patient Tubing*	3-551	10	25	
Tube Holder (Pack of 10)	3-552	1	2	
Fish Elbow (Pack of 10)*	3-530	2	5	
Bottle^	3-549	1	1	
Filter^	3-541	1	1	

Suggested quantities are based on the estimated number of patients to be treated per year and may need to be adjusted for individual hospitals.


\*Items labeled as single-use only  
 \*\*Hats may need to be changed every couple days to ensure they remain secure on the head  
 ^Replacement parts recommended yearly per Pumani CPAP device

## To Order Patient Accessories:

To request and order form and pricing information please contact us at:

 [sales@pumani.com](mailto:sales@pumani.com)

 [www.pumani.com](http://www.pumani.com)

 +1 (628) 236 7335

# Pumani

bubble CPAP



↑  
DO NOT CLOSE



# Pumani

bubble CPAP

Product Catalogue © 2018

For ordering information, please contact  
Hadleigh Health Technologies:



[www.pumani.com](http://www.pumani.com)



[info@pumani.com](mailto:info@pumani.com)



Hadleigh Health Technologies, LLC  
30 Castro Ave,  
San Rafael, CA 94901, USA  
+1 (628) 236 7335

